



**FlexAir Powered Air-Purifying Respirator (PAPR) Storage Battery Considerations**  
**User Notice**

ADG is posting this notice as a reminder to all users of, and entities that may store, the FlexAir PAPR System of battery maintenance requirements communicated in the *FlexAir Powered Air Purifying Respirator User Instruction* (P/N C-85090-00) provided with each unit.

The instruction, section "6.3 Storing the Respirator", informs:

*"The On-board rechargeable Li-ION Battery may require periodic charging during storage. For long term storage it is best to store the blower with batteries at a low charge level (~30% capacity). The battery should be charged about 10% every 6 months if not used for long periods. Follow the instructions included in this manual and those supplied with the AirBoss Defense Group Battery Charger."*

Long term storage, 6 months or greater, requires batteries be charged for forty-five minutes every six months. If battery maintenance has not been, or is not, performed, the battery's life may be significantly reduced.

Please contact Virginia Griffith at [Virginia.Griffith@ADG.com](mailto:Virginia.Griffith@ADG.com), or call 301-352-8800, with any questions.